

# Application Data: Chemical Industry

## Computrac® MAX® 4000XL Moisture/Solids Analyzer

MATERIAL	RESULT	COMPUTRAC®	STATISTIC	OVEN
Litter, Clay	% Moisture	9.5524	Mean	9.6043
		0.0592	Standard Deviation	0.0224
		6-7 minutes	Test Time	2 hours
Litter, Wood Granules	% Moisture	12.6464	Mean	12.9397
		0.0795	Standard Deviation	0.0708
		3 minutes	Test Time	2 hours
Powder, Epoxy	% Moisture	0.3034	Mean	0.3023
		0.0174	Standard Deviation	0.0267
		7 minutes	Test Time	2 hours
Powder, Iron	% Moisture	0.1324	Mean	0.1321
		0.0069	Standard Deviation	0.0074
		5 minutes	Test Time	2 hours



**ARIZONA INSTRUMENT LLC**

3375 N. Delaware Street, Chandler, AZ 85225 • (602) 470-1414 • (800) 528-7411 • [www.azic.com](http://www.azic.com)

*Exclusive Manufacturer of Computrac® Moisture Analyzers and Jerome® Toxic Gas Analyzers*