

# Application Data: Fly Ash

## Computrac<sup>®</sup> MAX<sup>®</sup> 5000 Moisture/Solids/Ash Analyzer

MATERIAL	RESULT	COMPUTRAC <sup>®</sup>	STATISTIC	Muffle Furnace ASTM C 311
Fly Ash Class F	% Moisture	0.1439	<b>Mean</b>	0.1180
		0.0073	<b>Standard Deviation</b>	0.0252
		2 minutes	<b>Test Time</b>	60 minutes (TCW)
Fly Ash Class F	% LOI	0.1423	<b>Mean</b>	0.1453
		0.0044	<b>Standard Deviation</b>	0.0124
		12-20 minutes	<b>Test Time</b>	60 minutes (TCW)



**ARIZONA INSTRUMENT LLC**

3375 N. Delaware Street, Chandler, AZ 85225 • (602) 470-1414 • (800) 528-7411 • [www.azic.com](http://www.azic.com)

*Exclusive Manufacturer of Computrac<sup>®</sup> Moisture Analyzers and Jerome<sup>®</sup> Toxic Gas Analyzers*