



COMPUTRAC®

by Arizona Instrument LLC

Computrac® MAX® 4000XL Highlights



A few unique features to the Computrac® MAX® 4000XL:

- MAX® 4000XL has a fan that comes on after the test is over to help cool the heating element down quicker. This produces a 25% quicker throughput of product.
- MAX® 4000XL is capable of allowing a customer to control the ramp rate. This comes in handy if a customer would like to know exactly when some of the other components in their product evaporate at temperature heating.
- MAX® 4000XL has the ability for Ethernet hook up for remote access through the internet. The software program is a little different and allows you to develop parameters for product.
- MAX® 4000XL allows the user to write the final test result formula. The Equation Editor calculates and displays results related to the parameters used to calculate moisture %.
- MAX® 4000XL stores 1000 of the last tests with 100 of the last graphs.
- MAX® 4000XL has an optional Bar Code Reader. If you are interested in identifying your sample using a Bar Code System, this optional feature is available.
- MAX® 4000XL has a Parameter Expert and Parameter Optimizer program that allows in house method development.

For additional information please contact your [Customer Service](#) or [Sales](#) Representative or call us at 800.528.7411.

[Brochure](#) | [Manual](#) | [Quote](#) | [Computrac® Tradeshows](#)



Updated 14 Oct 2011