




COMPUTRAC®

by Arizona Instrument LLC

Y990-0232 Web Server

The Computrac® Web Server Software is intended for use with the [Computrac® MAX® 4000XL](#) moisture | solids analyzer and will help access the instrument through an Intranet connection. The Web Server will allow remote viewing of the MAX® 4000XL's status, calibration information, test parameters, test results, and audit log from a networked PC. Once downloaded, data can be viewed and analyzed in Microsoft Excel, even after the instrument has been disconnected from the network.



The Computrac® Web Server is an optional feature available on MAX® 4000XL moisture | solids analyzers with firmware versions 1.0 and later. Contact Arizona Instrument today by phone at 800.528.7411 or 602.470.1414; or see the [User Manual](#) to learn more about this feature-packed, time-saving program.

System Requirements

Computrac® MAX® 4000XL:
Personal Computer:
Ports:

installed with firmware version 1.0 or later
installed with any Windows® based system with Internet Explorer 6.0 or later
one unused Ethernet port with assigned IP address