



COMPUTRAC®

by Arizona Instrument LLC

Y990-0214 Nitrogen Control Module Kit



The Nitrogen Control Module is an optional accessory that is used in conjunction with the [Nitrogen Purge Kit \(990-0046\)](#) to conserve nitrogen. The module turns on nitrogen to the instrument prior to a test and turns off the nitrogen after a test. This is especially helpful when using a limited supply of nitrogen (bottled nitrogen). The module also monitors the nitrogen pressure and alerts the user when the nitrogen source is depleted. The incoming pressure must be between 4 and 20 psi inclusive.

The Nitrogen Control Module is an optional feature available on [MAX® 4000XL](#) moisture and solids analyzers with firmware versions 2.1.01 and later. Contact Arizona Instrument today by phone at (800) 290-1414 or (602) 470-1414; or see the [User Manual](#) to learn more about this feature.

