



**COMPUTRAC®**

by Arizona Instrument LLC

## 990-0046 Nitrogen Purge Kit, MAX



In addition to the Nitrogen Purge Kit, the following materials are required for the nitrogen purge setup:

- Nitrogen source/cylinder or other supplied inert gas source
- Two-stage regulator capable of regulating pressure to 5 psi:

**For cylinder applications:** Multi-stage Gas Regulator for Compressed Air Cylinder, VWR #55850-420 or equivalent

**For supplied gas sources:** Matheson Gas Products Model 3471 Single Stage Line Regulator or equivalent

- Tubing with a 1/8" (.31 cm) inside diameter (I.D.) and, 1/4" (.625 cm) outside diameter (O.D.) - Vinyl tubing (VWR #63013-029), or Nalgene clear plastic tubing (Nalgene #8000-0020, VWR #63015-368)
- Connector, 1/8" I.D., 1/4" O.D.

Install the Nitrogen Purge Kit: [See User Manual](#)

This Kit (990-0046) May Be Used With: [Nitrogen Control Module Kit, Y990-0214](#)