



**COMPUTRAC®**  
by Arizona Instrument LLC

## COMPUTRAC® VAPOR PRO® Dry Air Generator Kit

The dry air generator kit is used in place of an on-site nitrogen source for Computrac® Vapor Pro® series moisture analyzers or for portable applications.



<b>Part Number</b>	Y990-0143 (110V)
<b>Dimensions</b>	6.00"L x 3.33"W x 6.00"H
<b>Weight</b>	3.00lbs
<b>Certifications</b>	UL, CSA
<b>Components</b>	Generator Drierite Canister Tygon Tubing Power Cord
<b>Options</b>	220V Converter (Part Number 990-0159)