

Wood Pellet Sample Steps 1-4

Step 1: Strainer, Mortar and Pestle, Wood Pellets



Step 2: Mash wood pellets



Step 3: Pour mashed pellets through strainer onto sample pan



Step 4: Stir to help push through mash pellets; ready for testing



ARIZONA INSTRUMENT LLC

3375 N. Delaware Street, Chandler, AZ 85225 • (602) 470-1414 • (800) 528-7411 • www.azic.com

Exclusive Manufacturer of Computrac® Moisture Analyzers and Jerome® Toxic Gas Analyzers