



Computrac® Printer Compatibility Vapor Pro® Series

Arizona Instrument LLC has verified that the following printers are compatible with Computrac® Vapor Pro® series moisture analyzers:

Firmware Revisions 1.12 and higher:

Printer Make and Model	Vapor Pro® Printer Driver
Canon BJC-250	Canon
Canon BJC-1000	Canon
Epson 440	Epson ESC/P2
Epson 880	Epson ESC/P2, Epson Raster, Epson C84
Epson C82	Epson ESC/P2, Epson Raster, Epson C84
Epson C84	Epson ESC/P2, Epson Raster, Epson C84
Epson C86	Epson ESC/P2, Epson Raster, Epson C84
Epson C88, C88+	Epson ESC/P2, Epson Raster, Epson C84
Seiko DPU-414	Seiko DPU

Firmware Revisions through 1.11A:

Printer Make and Model	Vapor Pro® Printer Driver
Epson Stylus C80, C82	Epson Old, Epson New
Epson 440, 740, 740i	IBM, Epson Old
Epson 880	Epson New
Canon BJC-250	Epson Old, Canon
Canon BJC-1000 and 2000	Canon
Seiko DPU (serial or parallel mode)	SII (for parallel interface)
IBM compatible ASCII/DOS	IBM
Epson compatible ASCII/DOS	Epson Old

Note: *With international installations, the end-user may see different characters on printed reports when compared to the characters seen on-screen. These differences are often caused by non-English character sets selected for printer output.*